

Programmation-connexions en un coup d'œil

| | NoahLink Wireless | ConnexxAir | Flex Strip Adaptateur | Adaptateur de Progr. 10 | Adaptateur de Progr. 312 | Adaptateur de Progr. 13 |
|-------------------|-------------------|------------|-----------------------|-------------------------|--------------------------|-------------------------|
| Insio Nx IIC | - | ✓ | - | - | - | - |
| Insio Nx CIC | - | ✓ | ✓ | - | - | - |
| Silk Nx | - | ✓ | ✓ | - | - | - |
| CROS Silk Nx | - | ✓ | ✓ | - | - | - |
| Pure 10 Nx | - | - | - | ✓ | - | - |
| Pure 312 Nx | ✓ | - | - | - | ✓ | - |
| Pure 13 Nx | ✓ | - | - | - | - | ✓ |
| Pure Charge&Go Nx | ✓ | - | - | - | - | - |
| Motion 13 Nx | ✓ | - | - | - | - | ✓ |
| Motion 13 P Nx | ✓ | - | - | - | - | ✓ |